**1**

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
**3**

- 1) Stick the tape (half width) in the deltoid muscle of the mid-humerus after shoulder joint.
- 2) Another tape (half width) stick from the same place over the front of delt.svalu towards the collarbone.
- 3) Another band (1.2 width) stick again from the same spot over the rear of the blade towards delt.svalu.

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**2**

- 1) Cut a piece of ribbon into a V-shaped stick from the head down along the cervical spine vertebra to the most prominent vertebra.
- 2) Next, stick the tape horizontally at the level of the lower neck,

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

**2**

- 1) Stick the tape from the wrist to the outside of the elbow.
- 2) Next stick the tape perpendicular to the first arm around.

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**2**

- 1) Stick a longer piece of tape from the side of the palm near the little finger to the bottom of the elbow.
- 2) Another tape stick around forearm.

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**2**

- 1) Stick the tape horizontally in strenuous back.
- 2) Next, stick the tape in X shape at the crossing point of the spine with the original tape.

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**EN**

**Please note:**

Apply to clean, dry and uninjured skin at least 1 hour before physical exercise. Do not use longer than 46 days after the application. When removing the tape (in the direction of hair), it is recommended to hold the skin with fingers to prevent the stretching. If you notice skin damage, remove the tape immediately. Keep out of reach of children! Do not expose the product to sunlight and high temperature.

**Recommendation:**

Before first use, contact your physiotherapist. Before applying, please read the attached instructions. Pregnant women should consult their doctor before using the tape. When applying, do not stretch the tape. Before use, stretch your muscles.

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